



# HYGIA

## Product Brochure Hygia BASE



Check In and Check Out  
To help with faster Contact Tracing



# HYGIA BASE - FEATURES



Completely  
Touchless Device



Light and Sound  
Notification



Plug and Play  
Installation



0.1 Second  
Response Time



Celsius and  
Fahrenheit Switch



USB and Battery  
Powered

# HYGIA BASE – SPECIFICATIONS

## SPECIFICATIONS

Accuracy	± 0.2 °C (34-45 °C, place it in the operating environment for 30 minutes before use)
Automatic Alarm	Flashing + “Buzzer” sound
Measurement	Measuring distance: 5cm - 10cm
Screen	Digital display
Charging Method	USB charging or battery (18650 Li-on)
Installation Method	Nail hook, Double-sided adhesive sticking, Bracket fixing, or Screw to Tripod
Environment Temperature	10°C - 40°C (Recommended 15°C - 35°C)
Infrared Measuring Range	0°C - 50°C
Response Time	0.5 seconds
Input	5V DC
Weight	380 gram
Dimensions	170 x 115 x 140 mm
Stand By	72 hours
Counting Capacity	999

**CAUTION:**


Devices are not be used under direct sunlight

# HYGIA BASE – OTHER DETAILS


## INSTALLATION METHODS



Wall Hook



Adhesive Tape



Tripod Stand

## COMPATIBILITY & CERTIFICATIONS



## WARRANTY & OTHER DETAILS

- All devices comes with **1-year** warranty on product and **6 months** warranty on accessories (tripod stand, batteries etc.).
- Compatible with SafeEntry.
- Lead Time: 2 weeks.
- Transportations cost are extra.

### CAUTION:

Devices are not be used under direct sunlight



# #HYGIAAGAINSTCOVID19

## ADDRESS

### SINGAPORE OFFICE

Datakrew Pte. Ltd. #02-26, 1 Cleantech  
Loop, Singapore 637141

### INDIA OFFICE

Datakrew Pvt. Ltd. HAL 2<sup>nd</sup> Stage,  
Indiranagar, Bangalore 560008

CAUTION:  
Devices are not be used under direct sunlight

CONFIDENTIAL DOCUMENT